

Economy Steel



Colour Options



Features

- > Removable Armrests
- > Fixed Wheels
- > Self propel and transit models
- > Robust steel frame

Product Code
SSP18E
STS18E

Superior Steel



Colour Options



Features

- > Half folding back
- > Luxury feel upholstery
- > Quick release rear wheels
- > Choice of metallic colours

Product Code
SUSP18
SUTS18

Lightweight Aluminium



Colour Options



Features

- > Our lightest self propel wheelchair
- > Outstanding value
- > Quick release rear wheels
- > LAW002 includes cable brakes

Product Code
LAWC001
LAWC002

Super Deluxe



Colour Options



Features

- > The ultimate in style, quality and durability
- > Adjustable for comfort
- > Several additional features

Product Code
XSSDSP
16/18/20/22
XSSDTS
16/18/20/22

Standard Steel



Colour Options



Features

- > Quick release rear wheels
- > Half folding back
- > Quick release footrests and armrests
- > Robust steel frame

Product Code
SSP18
STS18

Travel Chair 2



Colour Options



Features

- > Extremely lightweight and portable
- > Ideal for occasional trips
- > Comes with carry/storage bag
- > Quick release footrests

Product Code
TC002

XS Aluminium



Colour Options



Features

- > Anti-tip wheels with inbuilt stepper tube
- > Removable swing in/out footrests
- > Includes deluxe seat and back cushion
- > XSAWCTS18 includes cable brakes

Product Code
XSAWCSP18
XSAWCTS18

SPIRIT



Colour Options



Features

- > Stylish and lightweight frame design
- > One touch armrests
- > Puncture proof tyres
- > Long reach brake handles

Product Code
XSES18

PowerSTROLL



Features

- > Maximum speed 4mph
- > Maximum range 10 miles
- > Variable speed dial and thumb throttle
- > Easy to install and remove
- > 3 models available including single and dual wheel

Product Code

- PWCPP004
- PWCPP005
- PWCPP006

> Accessories

- > Wheelchair Bags and Holdalls
- > Wheelchair Clothing
- > Wheelchair Cushions
- > Telescopic Ramps
- > Rollup Ramps
- > Lap Belts
- > Wheels and Castors
- > Handrims and Tyres
- > Wheelchair Umbrella
- > Puncture Proof Tyres



Code	Dimensions cm (L x W x H)			Total Weight kg (lb)	Carry Format Weight kg (lb)	Weight Capacity kg (st)
------	------------------------------	--	--	----------------------------	-----------------------------------	-------------------------------

Steel Wheelchairs						
SSP16	104	61	95	17.8 (39.1)	10.2 (22.4)	115 (18)
SSP18	104	66	95	18.3 (40.2)	10.5 (23.1)	115 (18)
SSP18E	104	69	95	19 (41.8)	15 (33)	115 (18)
SUSP18	104	66	95	18.3 (40.2)	10.5 (23.1)	115 (18)
SSP22	104	76	95	20.1 (44.2)	12.3 (27)	130 (20)
STS16	97	58	95	15.6 (34.3)	10 (22)	115 (18)
STS18	97	63	95	15.1 (33.2)	10.5 (23.1)	115 (18)
STS18E	97	63	95	17 (37.4)	13 (33)	115 (18)
SUTS18	97	63	95	15.1 (33.2)	10.5 (23.1)	115 (18)
STS22	97	73	95	17.9 (39.5)	12.3 (27.1)	130 (20)

Aluminium Wheelchairs						
TC002	87	56	90	8.6 (19)	7.4 (16.3)	115 (18)
LAWC001	106	61	90	14.5 (31.9)	8.5 (18.7)	115 (18)
LAWC002	96	55	91	12.4 (27.3)	8.5 (18.7)	115 (18)
XSAWCSP18	103	63	93	16.8 (37)	8.3 (18)	115 (18)
XSAWCTS18	92	63	93	14.7 (32)	8.3 (18)	115 (18)
XSSDSP16	105	61	100	18.2 (40.1)	9.3 (20.5)	135 (21)
XSSDSP18	105	66	100	18.4 (40.6)	9.5 (21)	135 (21)
XSSDSP20	105	71	100	18.6 (41)	9.7 (21.3)	135 (21)
XSSDSP22	105	76	100	18.8 (41.4)	9.9 (21.8)	135 (21)
XSSDTS16	102	57	100	16.9 (37.3)	10.2 (22.5)	135 (21)
XSSDTS18	102	62	100	17.1 (37.7)	10.4 (22.9)	135 (21)
XSSDTS20	102	67	100	17.3 (38.1)	10.6 (23.4)	135 (21)
XSSDTS22	102	72	100	17.5 (38.6)	10.8 (23.8)	135 (21)
XSES18	97	71	91	14.9 (33)	9.2 (20.3)	115 (18)

Powerstroll					
PWCPP004	N/A		12.4 (27)	5.5 (12)	135 (21)
PWCPP005	N/A		20 (44)	10.8 (24)	135 (21)
PWCPP006	N/A		21 (46)	11.8 (26)	135 (21)

drive

MEDICAL DESIGN & MANUFACTURING



“ Our Enigma range of wheelchairs guarantee a focus on Comfort, Style, Performance and Reliability ”

The Enigma Wheelchair range



We understand the true value of innovation within this market place; and know that our Wheelchair development will have a direct influence on improving an individual's life. For this reason ALL the Wheelchairs in our range share the following features: -

- > Foldable for storage
- > Removable footrests
- > Puncture proof front castors
- > Height adjustable footrests
- > Comfortable Upholstery
- > Tool kit included
- > Simple to operate
- > Quick release components
- > Parking brakes

Leaflet Code: LL089



ENIGMA

The best Wheelchair Range ever from



www.drivemedical.co.uk

All sizes and weight stated are nominal. All errors and omissions excepted. Due to the limitations of photography and the printing process colours shown may not be 100% accurate.