

Libre LX

UNBEATABLE CONVENIENCE
AND VALUE.

LIBRE LX

The **Libre LX's** attractive, compact design makes it convenient and easy to take along with you. Standard features like a removable battery pack, easily adjustable tiller, and off-board battery charger make the Libre LX fun and simple to use.



BUILT WITH PRIDE

Pride[®]
Mobility Products Ltd.

(Tel) 01869 324 600 • (Fax) 01869 323 070 • www.pridemobility.com

Libre LX

Weight capacity:	18 Stone (113 kg)
Maximum speed:	Up to 6 km/h (4 mph)
Ground clearance:	4 cm (1.5")
Turning radius:	110.5 cm (43.5")
Length:	108 cm (42.625")
Width:	49.5 cm (19.5")
Front tyres:	5 cm x 18 cm (2" x 7")
Rear tyres:	6 cm x 20 cm (2.5" x 8")
Seat dimensions:	43 cm x 41 cm (17" x 16")
Anti-tips:	Rear
Standard seating:	Compact foldable
Battery charger:	O -board, 2 amp
Battery requirements:	(2) 12V deep cycle
Battery size:	12 AH
Battery pack weight:	7.5 kg (16.5 lbs.)
Seat weight:	10 kg (22 lbs.)
Heaviest piece when disassembled:	23 kg (50.5 lbs.) (base weight without seat and batteries)
Total overall weight:	42 kg (92.5 lbs.)

Color:



Red



Blue



The information contained herein is correct at the time of publication; due to the manufacturers commitment to constant improvement and development, they reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL
(Tel) 01869 324 600 • (Fax) 01869 323 070 • www.pridemobility.com

Pride[®]
Mobility Products Ltd.