



# RANGER



The Pride® Ranger provides optimal power, superior performance, unmatched reliability and superb style. With its powerful dual motors, the Ranger is the ultimate all-terrain personal mobility vehicle.

*Stylish design  
&  
premium  
performance  
at an excellent  
value!*

**RANGER**



Tel: 01869-324600  
Fax: 01869-323070  
[www.pridemobility.com](http://www.pridemobility.com)



**Pride**  
Mobility Products Ltd.

## Specifications

Weight capacity	181 kg (28.5 STONE)
Maximum speed	13 km/h (8 mph)
Ground clearance	18 cm (7")
Turning radius	196 cm (77")
Length	161 cm (63.5")
Width	79 cm (31")

## Tyres

Front	15 cm x 33 cm (6" x 13") pneumatic
Rear	15 cm x 33 cm (6" x 13") pneumatic
Anti-tips	N/A
Standard seat	Highback Grey vinyl
Power elevating seat	Option: yes
Suspension	Rear
Battery charger	Off-board, 8 amp
Per-charge range	Up to 40 km (25 miles)
Battery requirements	(2) 12V deep-cycle
Battery Size(s)	100AH (Group 27)

## Weight

Battery	29 kg (63.5 lbs) each
Total weight without batteries	98 kg (216 lbs)
Heaviest piece when disassembled	No frame disassembly

## Colours



Candy  
Apple Red



Black  
Onyx

*"The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice".*

*"Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement".*

## Features

- 200 amp Curtis microprocessor based PMV controller system provides optimal levels of torque and power
- Two rugged independent in-line gear motors provide positive traction
- Variable speed up to 13 km/h (8mph)
- Heavy duty 15 cm x 33 cm (6" x 13") all-terrain tyres with Pride "Micro-Mag" alloy wheel assemblies for superb, stable outdoor performance
- Seating features user adjustable sliding mechanism

## Options

- Lap belt
- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Walker holder
- Cane clip attachment
- Golf bag attachment

